1	Intended use
2	Summary and explanation
3	Principle of the assay
4	Reagents
5	Caution
6	Storage and stability
7	Specimen collection and preparation
8	Methodology
9	Ring measurement and result processing
10	Limitations of procedure
11	Expected values
12	Performance characteristics
13	Bibliography
14	Summary of procedure

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For *in-vitro* diagnostic use only

Product Code: RC001.1/2/3

Product manufactured by:
The Binding Site Group Ltd, PO Box 11712, Birmingham B14 4ZB, U.K.
www.bindingsite.co.uk
Telephone: +44(0)121 436 1000
Fax: +44(0)121 430 7061
e-mail: info@bindingsite.co.uk

FDA (USA) Information
Analyte ID 1046
Test System 61450
Complexity Cat. High

1 INTENDED USE

This kit is intended for measuring human haemolytic complement activity by the classical pathway in human serum, as an aid in diagnosis of immunological disorders, especially those associated with complement component deficiencies.

2 SUMMARY AND EXPLANATION

The classical complement pathway is a cascade of reactions involving 11 proteins, which leads to lysis of cell membranes and ultimately to cell death.

Activation of the classical complement pathway is brought about by the interaction of antibody-antigen complexes with complement C1, which then activate C4 causing in turn activation of C2, C3 and C5. Activated C5 (C5b) binds C6, C7, C8 and 12 to 18 molecules of C9 to form a large complex known as the membrane attack complex (MAC). It is this complex which disrupts cell membranes leading to cell death. (refs. 1-6)

Decreased complement activity may be caused by deficiencies of any of the individual complement components, hereditary or acquired, and can be associated with hereditary complement deficiencies, glomerulonephritis, systemic lupus erythematosus (SLE), vasculitis and a number of non-immunological conditions. (ref. 7)

3 PRINCIPLE OF THE ASSAY

The method involves serum containing active complement components diffusing radially at 2-6°C from a cylindrical well through an agarose gel containing sheep erythrocytes sensitised with a rabbit anti-sheep erythrocyte antibody, termed 'Haemolysin'. During this diffusion period complement C1 is fixed by the Haemolysin. When incubated at 37°C the complement cascade proceeds, leading to lysis of the sheep erythrocytes, resulting in a clear zone, the size of which is dependent on the original complement activity in the serum.

By measuring the zones of lysis produced by a number of sera of known complement activity a calibration curve may be constructed by plotting the diameter against complement activity on semi-log graph paper. The complement activity of an unknown sample may then be determined by measuring the zone of lysis and reading off the calibration curve providing the ring quality is acceptable (see limitations of procedure, section 10).

4 REAGENTS

- 4.1 Haemolytic complement plates (supplied in resealable foil pouches). These contain sensitised sheep erythrocytes. Up to fourteen samples can be run per plate (including calibrators and controls). Preservative: 0.02% sodium azide.
- 4.2 Calibrators. These are supplied lyophilised. The complement activity is given on the vial label.
- 4.3 Controls. These are supplied lyophilised. The expected complement activity is given on the vial label.

5 CAUTION

All donors of human materials supplied in this kit have been serum tested and found negative for Hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus and human immunodeficiency virus (HIV1 and HIV2). However, these tests cannot guarantee the absence of infective agents. Proper handling and disposal methods should be established as for all potentially infective material and only personnel adequately trained in such methods should be permitted to perform the procedures.

The plates contain sodium azide as a preservative and must be handled with caution-do not ingest or allow contact with the skin or mucous membranes. If contact does occur wash with a large volume of water and seek medical advice. Explosive metal azides may be formed with the lead and copper plumbing; on disposal of reagent, flush with a large volume of water to prevent azide build up.

Adherence to the given procedure is recommended. The validity of results obtained using methods other than those stated cannot be guaranteed. Reagents from different batch numbers of kits are NOT interchangeable. If large numbers of tests are performed, care should be taken to ensure that all reagents are from the same batch.

6 STORAGE AND STABILITY

The unopened kits should be stored at 2-8°C and can be used until the expiry date given on the kit box label. DO NOT FREEZE. The haemolytic complement plates have a shelf life limited by the survival of the sheep erythrocytes and should be kept refrigerated wherever possible. The agarose gel will be destroyed by freezing, therefore keep plates away from cooling elements in refrigerators. Unopened, plates should be stored flat and upside down (pouch label uppermost) to prevent condensation accumulating in the wells. Handle plates with care to prevent gel damage.

Once reconstituted, calibrators and controls should be used immediately. If storage of reconstituted material is required, e.g. for use with part-used plates, aliquots should be frozen, at -20°C or below.

7 SPECIMEN COLLECTION AND PREPARATION

Complement activity may be activated in improperly handled samples leading to erroneous results. The correct collection, processing, storage and shipment of specimens is therefore essential.

Use fresh serum samples. Microbially contaminated, haemolysed and very lipaemic serum samples or those containing particulate matter should not be used as incomplete lysis may occur. Blood samples should be collected by venepuncture, and allowed to clot at room temperature (18-24°C) for approximately 60 minutes. Centrifuge the sample and transfer to a clean tube. Store for a maximum of 4 hours at 2-8°C or on ice prior to assay. For longer term storage freeze samples at -70°C or below. Frozen specimens should be thawed slowly on ice and tested as soon as possible after thawing. Repeated freezing and thawing should be avoided.

8 METHODOLOGY

(A summary of the entire procedure is given at the end of this instruction leaflet.)

8.1 Materials provided

- 8.1.1 1 x Total Haemolytic Complement Plate
in foil pouch (2 plates in kit RC001.2, 3 plates in kit RC001.3)
- 8.1.2 8 x Gel Dividers
- 8.1.3 1 x Total Haemolytic Complement Calibrator
(2 in 2 plate kit, 3 in 3 plate kit)
- 8.1.4 1 x Total Haemolytic Complement Control Serum
(2 in 2 plate kit, 3 in 3 plate kit)
- 8.1.5 1 x Instruction leaflet
- 8.1.6 Semi-log graph paper

8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples, e.g. samples tubes, centrifuge etc.
- 8.2.2 Pipettes for accurate reconstitution of calibrators and controls and for dilution of calibrators.
- 8.2.3 Micropipettes for sample application. These should be capable of accurately delivering 5µL volumes. Binding Site Micropipettes (code AD041) or 'Hamilton' syringes are recommended.
- 8.2.4 Jewellers' Eyepiece (code AD040) or Digital RID Plate Reader (AD400) for magnifying and accurately measuring the zones of lysis to 0.1mm.
- 8.2.5 Physiological saline.
- 8.2.6 Distilled water.

8.3 Reagent preparation

8.3.1 RID plate

To avoid contamination of the gel, plates should be used in a dust-free environment. Take the plate from the foil pouch and remove the lid. If condensation is visible the plate should be kept upside down until the lid has been removed to prevent droplets falling onto the gel. Check the plate to ensure that no damage has occurred in storage or transit, e.g. splits in the gel. Leave the plate open for 10-15 minutes (or longer if necessary) at room temperature to allow any condensation in the wells or on the gel surface to evaporate. Samples should never be applied to wells in which moisture is still visible.

Plate partitioning: The plates may be partitioned into up to four sections using the gel dividers provided. Each divider should be positioned carefully on the gel, cutting edge downward, with the stabilising arm resting on the central plate label. Press firmly on the arm to cut the gel and leave in position.

Plate partitioning is recommended if the plate is to be used on two separate occasions. Plates must be partitioned prior to use. After initial use, partitioned plates should be resealed in their foil pouches and stored at 2-8°C with the gel divider(s) in place. Store partitioned plates right side up.

8.3.2 Calibrator

The lyophilised calibrators should be reconstituted with distilled water pre-cooled to 2-8°C using the volume indicated on the vial labels. Before use, all material in the bottle, including any adhering to the bung must be completely dissolved by inversion. Prepare a 1/2 dilution by mixing 100µL of calibrator with 100µL of physiological saline which has been pre-cooled to 2-8°C and a 1/4 dilution by mixing 100µL of calibrator and 300µL of pre-cooled saline. The calibrator and its dilutions must be kept as cool as possible during preparation to maintain the integrity of the complement components.

8.3.3 Control

The lyophilised control should be reconstituted with distilled water pre-cooled to 2-8°C using the volume indicated on the vial label. It should be mixed gently by inversion until the contents are completely dissolved.

8.3.4 Sample

Samples should be applied neat. However, if a dilution is required use physiological saline, which has been pre-cooled to 2-8°C.

8.4 Application of calibrators and samples

For best results it is recommended that all reagents, samples and plates etc, be kept cool throughout the application procedure. The use of a cold block will enable the samples to be applied to pre-cooled plates.

The calibrators, control and test samples should be gently mixed immediately before use. 5µL of each of the three calibrator dilutions should be pipetted into three separate wells. The remaining wells should then be filled with 5µL aliquots of samples and control as required. Plates should not be left open for long periods during calibrator / sample application, as this will cause excessive drying of the gel.

8.5 Diffusion

These haemolytic complement plates can be used to perform either a same day assay or an overnight assay. For the same day assay, the samples should be allowed to diffuse at 2-6°C for at least 6 hours to ensure the formation of sufficiently large lytic zone diameters (see limitation 10.2).

When being used to perform an overnight assay, diffusion can be carried out for up to 18 hours. If diffusion is allowed to continue for longer than 18 hours multiple rings and / or incomplete lysis may occur. For both assays the diffusion step should be carried out at 2-6°C but for overnight assay it is especially important that the temperature be kept as close to 4°C as possible to prevent degradation of the more fragile complement components. This will ensure that the final zones of lysis are sharp with well-defined boundaries. **Incubation temperatures >6°C during the diffusion period may prevent complete lysis.**

8.6 Incubation

Following the required period of diffusion, plates should be incubated at 37°C for a minimum of 30 minutes. It is important that the gels do not dry out and to minimise this plates should be resealed in their foil pouches.

8.7 Quality control


The control should, following reconstitution, be treated exactly like a test sample. Values for the control should be within 20% of the value stated on the vial label.

9 RING MEASUREMENT AND RESULT PROCESSING

After incubation the diameters of the zones of lysis should be measured to the nearest 0.1mm using a jewellers' eyepiece or RID reader. When reading with an eyepiece, use bright side lighting and a dark background. Using the semi-log graph paper provided, plot the diameters (in mm) of the rings formed by the three calibrators along the linear y-axis against their % complement activity (100% for the neat calibrator and 50% and 25% for the two dilutions) along the log x-axis. The three points should then be joined together as shown below. From this curve, % complement activity values for the test samples and controls can be read, and these then converted to CH100 units/mL values.

Example: The following calibration curve was obtained on a plate where diffusion was allowed to continue for 18 hours at 2-6°C. The value of the neat calibrator (100% activity) was 1171 CH100 units/mL.

A test sample applied neat as recommended gave a ring diameter of 8.0mm. From the calibration curve this corresponds to an activity of 45% and hence a complement activity of 527 CH100 units/mL (1171 x 0.45).

10 LIMITATIONS OF PROCEDURE

- 10.1 The smaller rings obtained by the 6hr assay will give a "flatter" graph making the results less accurate than those obtained with the overnight method.
- 10.2 For samples where there is deficiency in overall complement activity or of a single complement component, the resulting zones of lysis may contain unlysed cells or show multiple rings. Where partial lysis or poor ring quality is observed the quantitative value obtained may be inaccurate. In such cases, a repeat sample should be requested. If the same result is obtained, then samples should be further investigated for dysfunction of the complement pathway.
- 10.3 Haemolysed samples or sera from rheumatoid arthritis patients may also show poor or multiple rings.
- 10.4 If two or more plates are to be used simultaneously a separate calibration curve should be constructed for each plate.
- 10.5 Reuse of part-used plates: The unused wells of part-used plates can be used for assaying further samples on a subsequent occasion providing the plate shows no signs of deterioration. A new calibration curve is required each time. The zones of lysis around previously used wells may increase in size due to further incubation at 37°C. This should not be considered as deterioration and will not affect the performance of the plate.
- 10.6 Calibrators and controls may give rings of poorer quality than fresh samples. This will not affect the accuracy of the assay.

10.7 Trouble shooting

Problem	Possible Cause	Suggested Action(s)
A. No ring for calibrators or control	i) Calibrator/control omitted.	Repeat assay.
	ii) Incubation temperature during diffusion period >6°C.	Repeat assay.
B. No ring for test sample	i) Sample omitted.	Repeat assay.
	ii) Incubation temperature during diffusion period >6°C.	Repeat assay.
	iii) Sample haemolysed during collection or storage.	Review sample collection and storage. Repeat with a fresh sample.
	iv) Haemolytic complement activity low due to classical pathway dysfunction or poor storage/handling of sample.	Review sample storage and handling. Repeat assay if necessary.
C. Non-circular rings	i) Poor sample application.	Repeat assay.
	ii) Gel dried out before use, due to storage damage or expired product.	Repeat assay using new plate.
	iii) Gel dried out during sample application or incubation.	Repeat assay using new plate, minimising time the plate is left open. Incubate with lid on tight and seal in foil pouch.
D. Clear zones on gel prior to use	Plate has been frozen, dried out, expired and/or been contaminated.	Review storage and repeat assay using new plate.
E. Incomplete lysis and/or multiple rings	i) Haemolytic complement activity low due to classical pathway dysfunction or factor deficiency.	Confirm with other clinical findings.
	ii) Sample haemolysed during collection or storage.	Review sample collection and storage. Repeat with a fresh sample.
	iii) Rheumatoid arthritis patient sera.	Check patient's clinical history.
	iv) Incubation temperature during diffusion period >6°C.	Repeat assay.
F. Weak, pitted gel	Plate has been frozen.	Repeat assay using new plate. Review storage.

- 10.8 If an unexpected result is obtained, the assay should be repeated, preferably with a fresh sample.
- 10.9 This kit is used to aid diagnosis only. An abnormal result may suggest certain diseases, which must be confirmed by clinical findings and other serological tests. The results obtained from this assay are not diagnostic proof of the presence or absence of disease.
- 10.10 For FDA (USA) information see front page.

If problem persists, please refer to supplier.

11 EXPECTED VALUES

Normal adult serum complement activities (determined using these kits):

	Median activity (CH100 units/mL)	95 Percentile range (CH100 units/mL)
Number of samples (n) = 38	728	392-1019

Decreased complement activity can be associated with various clinical conditions (see Section 2). Expected values are affected by age, sex, diet, geographical location, pregnancy, illness and other factors. To overcome differences in sample collection/preparation and assay procedure it is strongly recommended that each laboratory establish its own normal ranges. Samples with a complement activity below the 95 percentile range may be complement deficient and should be investigated further. Where partial lysis or poor ring quality is observed the quantitative value obtained may be inaccurate. In such cases, a repeat sample should be requested. If the same result is obtained, then samples should be further investigated for dysfunction of the complement pathway.

12 PERFORMANCE CHARACTERISTICS**12.1 Precision**

Intra-assay precision was determined by measuring three serum preparations using the 18 hour assay method on five THC plates from a typical batch. Inter-assay precision was determined by measuring a serum preparation using the 18 hour assay method on five THC plates from different batches.

INTRA-ASSAY PRECISION		
	Mean conc. (CH100 units/mL)	CV (%)
Sample 1	1050	0
Sample 2	840	3.9
Sample 3 (Control)	590	10.3

INTER-ASSAY PRECISION		
	Mean conc. (CH100 units/mL)	CV (%)
Sample A	627	1.6

12.2 Comparison

The relative specificity, sensitivity and agreement have been determined against a commercially available reference method using 79 clinical sera and 40 normal adult donor sera. Samples were regarded as negative on the Binding Site kit if their CH100 activity was below the 95 percentile range.

Binding Site assay		Competitor assay	
		+	-
	+	29	11
	-	1	78
Relative sensitivity		96.7%	
Relative specificity		87.6%	
Relative agreement		89.9%	

13 BIBLIOGRAPHY

- 13.1 Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
- 13.2 Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin. Immunol.* 5, 83-89.
- 13.3 Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
- 13.4 Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
- 13.5 Kolb WP et al. (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J Exp. Med.* 135, 549-566.
- 13.6 Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
- 13.7 Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

14 SUMMARY OF PROCEDURE

- 14.1 Reconstitute calibrator and control with distilled water precooled to 2-8°C.
- 14.2 Prepare 1/2 and 1/4 dilutions of the calibrator using cold saline solution.
- 14.3 Allow condensation to evaporate from the plate(s).
- 14.4 Apply the calibrator, calibrator dilutions, control and samples to the plate in 5µL volumes.
- 14.5 Replace lid and allow diffusion to occur for between 6 and 18 hours at 2-6°C.
- 14.6 Incubate at 37°C for a minimum of 30 minutes.
- 14.7 Measure zones of lysis.
- 14.8 Plot calibration curve and read off results.

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Bestell-Nr.: RC001.1/2/3

In England hergestellt von:

The Binding Site Group, P.O. Box 11712, Birmingham B14 4ZB, United Kingdom.
www.bindingsite.co.uk

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13 REFERENZEN

- Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
- Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin.Immunol.* 5,83-89.
- Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
- Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
- Kolb WP et al. (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J Exp. Med.* 135, 549-566.
- Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
- Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

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13 BIBLIOGRAPHIE

- 13.1 Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
 13.2 Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin. Immunol.* 5, 83-89.
 13.3 Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
 13.4 Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
 13.5 Kolb WP et al. (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J Exp. Med.* 135, 549-566.
 13.6 Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
 13.7 Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

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- Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
- Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin. Immunol.* 5,83-89.
- Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
- Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
- Kolb WP et al. (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J Exp. Med.* 135, 549-566.
- Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
- Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

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13 BIBLIOGRAFIA

- 13.1 Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
- 13.2 Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin. Immunol.* 5, 83-89.
- 13.3 Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
- 13.4 Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
- 13.5 Kolb WP *et al.* (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J. Exp. Med.* 135, 549-566.
- 13.6 Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
- 13.7 Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

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13 BIBLIOGRAFIA

- Janeway CA & Travers P (1994). Immunobiology: The immune system in health and disease. Blackwell Scientific Publications, Oxford.
- Tomlinson S (1993). Complement defence mechanisms. *Curr. Opin. Immunol.* 5, 83-89.
- Frank MM & Fries LF (1991). The role of complement in inflammation and phagocytosis. *Immunol. Today* 12, 322-326.
- Poddack ER & Tschopp J (1982). Polymerisation of the ninth component of complement (C9): the formation of poly-C9 with a tubular ultrastructure resembling the membrane attack complex of complement. *Proc. Natl. Acad. Sci.* 79, 574-578.
- Kolb WP et al. (1972). Molecular analysis of the membrane attack complex mechanisms of complement. *J. Exp. Med.* 135, 549-566.
- Muller-Eberhard HJ (1975). Complement. *Ann. Rev. Biochem.* 44, 697-724
- Morgan BP & Walport MJ (1991). Complement deficiency and disease. *Immunol. Today* 12, 301-306.

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