1	Intended use
2	Summary and explanation
3	Principle of the assay
4	Reagents
5	Caution
6	Storage and stability
7	Specimen collection and preparation
8	Methodology
9	Ring measurement and result processing
10	Limitations of procedure
11	Expected values
12	Performance Characteristics
13	Bibliography
14	Summary of procedure
15	RID reference table

Deutsch

Français

Español

HUMAN COMPLEMENT C1q BINDARID™ RADIAL IMMUNODIFFUSION KIT

For *in vitro* diagnostic use only

Product code: RN020.3

BINDARID™ is a trademark of The Binding Site Group Ltd., Birmingham, U.K.

Product manufactured by:

*The Binding Site Group Ltd., PO Box 11712, Birmingham B14 4ZB, U.K.
www.bindingsite.co.uk*

Telephone: +44 (0) 121 436 1000

Fax: +44 (0) 121 430 7061

e-mail: info@bindingsite.co.uk

FDA (USA) Information

Analyte ID Code: 1064

Test System ID Code: 61073

Complexity Cat: High

1 INTENDED USE

This kit is intended for measuring human C1q in serum as an aid in the diagnosis and treatment of systemic lupus erythematosus (SLE).

2 SUMMARY AND EXPLANATION

C1q is a 400kD hexameric gamma-2 protein that is a subunit of C1, the first component of complement. The binding of two or more of C1q's six globular domains initiates the classical pathway of complement activation. C1q binds readily to the CH2 domains of aggregated IgG molecules in an immune complex or the CH3 domains of a single IgM molecule whose conformation has been altered following antigen binding. It can also bind directly to certain micro-organisms and mycoplasmas. The multivalent binding of C1q is believed to lead to a conformational change in the C1q complex, activating C1r then C1s and thereby initiating the classical complement pathway. The collagen-like tail domain of C1q (which is only exposed once C1 inactivator dissociates C1r₂ and C1s₂ from the C1q-activator complex) increases the phagocytosis of particles by monocytes and macrophages. Serum levels of C1q are reduced in immune complex disease, SLE and meningitis. Hereditary deficiency is also known (refs 1-4).

Radial immunodiffusion (RID) is a technique that is routinely used for measuring the concentration of various soluble antigens in biological fluids. It is principally derived from the work of Fahey & McKelvey (ref. 5) and Mancini *et al* (refs. 6 & 7).

3 PRINCIPLE OF THE ASSAY

The method involves antigen diffusing radially from a cylindrical well through an agarose gel containing an appropriate mono-specific antibody. Antigen-antibody complexes are formed which, under the right condition, will form a precipitin ring. The ring size will increase until equilibrium is reached between the formation and breakdown of these complexes, this point being termed 'completion'. At this stage, a linear relationship exists between the square of the ring diameter and the antigen concentration. By measuring the ring diameters produced by a number of samples of known concentration, a calibration curve may be constructed. The concentration of the antigen in an unknown sample may then be determined by measuring the ring diameter produced by that sample and reading off the calibration curve.

There are three different procedures that may be used with this kit (see Section 8.4). Procedures ONE and TWO require that rings are measured at completion. A linear calibration curve is constructed for Procedure TWO, whereas for Procedure ONE a reference table (based upon the ideal linear calibration curve) is provided, which converts ring diameters directly to protein concentrations. Using Procedure THREE, ring diameters are measured before completion; the calibration curve produced will be non-linear.

4 REAGENTS

- 4.1 RID plates (supplied in foil pouches). These contain monospecific antibody to C1q in agarose gel. Up to fourteen samples can be run per plate (including calibrators). Preservatives: 0.099% sodium azide, 0.1% E-amino-n-caproic acid (EACA), 0.01% thiomersal (sodium ethylmercurithiosalicylate), 0.01% benzamidine.
- 4.2 Calibrators. These are supplied lyophilised as a set of three containing high, medium and low concentrations of C1q. The concentrations of C1q given on the vial labels have been obtained by comparison with a commercially available standard; however in the absence of international agreement the accuracy of the available standard cannot be guaranteed. Preservatives: 0.099% sodium azide, 0.1% EACA, 0.01% benzamidine.
- 4.3 7% Bovine Serum Albumin (BSA) solution. This is supplied in stabilised liquid form and is included for use as a diluent. Preservative: 0.099% sodium azide, 0.1% EACA, 0.01% benzamidine.
- 4.4 Control. This is supplied lyophilised. The expected C1q concentration is marked on the bottle label. Preservatives: 0.099% sodium azide, 0.1% EACA, 0.01% benzamidine.
- 4.5 Distilled water. For reconstituting the lyophilised calibrators and control. Preservative: 0.099% sodium azide.

5 CAUTION

All donors of human serum supplied in this kit have been serum tested and found negative for Hepatitis B surface antigen (HBsAg) and antibodies to human immunodeficiency (HIV1 and HIV2) and Hepatitis C virus. However, these tests cannot guarantee the absence of infective agents. Proper handling and disposal methods should be established as for all potentially infective material and only personnel adequately trained in such methods should be permitted to perform the procedures.

The plates, calibrators, controls and BSA contain 0.099% sodium azide and/or 0.01% thiomersal as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. If contact does occur wash with a large volume of water and seek medical advice. On disposal of reagent flush with a large volume of water, to prevent azide build up.

Adherence to the given procedure is recommended. The validity of results obtained using methods other than those stated cannot be guaranteed.

Reagents from different batch numbers of kits are NOT interchangeable. If large numbers of tests are performed, care should be taken to ensure that all reagents are from the same batch.

6 STORAGE AND STABILITY

The unopened kits should be stored at 2-8°C and can be used until the expiry date given on the kit box label. DO NOT FREEZE. The expiry dates of individual components are given on the component labels. RID plates should be stored at 2-8°C and are damaged by temperature extremes. Freezing will destroy the gel, therefore RID plates should be kept away from cooling elements in refrigerators. High temperatures should also be avoided as this will result in moisture loss from the gel, affecting performance. Unopened plates should be stored flat and upside down (pouch label uppermost) to prevent condensation accumulating in the wells. Handle plates with care to prevent gel damage.

Unopened calibrators and controls should be stored at 2-8°C. Once reconstituted they are stable for at least one week at 2-8°C, but for longer storage they should be aliquoted and frozen (-20°C or below). All other reagents should be stored at 2-8°C.

7 SPECIMEN COLLECTION AND PREPARATION

Use fresh or deep frozen (-20°C or below) serum samples. Microbially contaminated, haemolysed and very lipaemic serum and samples containing particulate matter should not be used. Blood samples should be collected by venepuncture, allowed to clot naturally and the serum separated as soon as possible to prevent haemolysis. The serum may be stored at 2-8°C for up to 48 hours prior to assay, or for prolonged storage, aliquoted and kept at -20°C or below. Repeated freezing and thawing should be avoided.

The BSA included in the kit should be used as diluent when required, as this will maintain the viscosity of the material. Results can therefore be accurately compared with the calibrator which has a similar viscosity to normal serum.

8 METHODOLOGY

(A summary of the entire procedure is given at the end of this instruction leaflet).

8.1 Contents:

- 8.1.1 3 x Human Complement C1q NL Bindarid (radial immunodiffusion plates in foil pouches)
- 8.1.2 8 x Gel Dividers
- 8.1.3 3 x Human C1q NL Calibrator (lyophilised)
- 8.1.4 1 x 5mL 7% BSA Solution
- 8.1.5 1 x Human C1q Control Serum (lyophilised)
- 8.1.6 1 x 5mL Distilled water
- 8.1.7 1 x Instruction leaflet, including RID reference table

8.2 Materials required but not provided:

- 8.2.1 Equipment for collection and preparation of test samples, e.g. sample tubes, centrifuge etc.
- 8.2.2 Pipettes for accurate reconstitution of calibrators and control and dilution of samples.
- 8.2.3 Micropipettes for sample application. These should be capable of accurately delivering 10µL volumes. Binding Site Micropipettes (code AD041) or 'Hamilton' syringes are recommended.
- 8.2.4 Jeweller's eyepiece (code AD040) or digital RID plate reader (code AD400) for magnifying and accurately measuring the precipitin ring diameters to 0.1mm.
- 8.2.5 Graph paper.

8.3 Reagent preparation

8.3.1 RID plate(s)

To avoid contamination of the gel, plates should be used in a dust-free environment. Take the plate from the foil pouch and remove the lid. If condensation is visible the plate should be kept upside down until the lid has been removed to prevent droplets falling onto the gel. Check the plate to ensure that no damage has occurred in storage or transit, e.g. splits in the gel. Leave the plate open for 10-15 minutes (or longer if necessary) at room temperature to allow any condensation in the wells or on the gel surface to evaporate. Samples should never be applied to wells in which moisture is still visible.

Plate partitioning: The plates may be partitioned into up to four sections using the gel dividers provided prior to use. Each divider should be positioned carefully on the gel, cutting edge downward, with the stabilising arm resting on the central plate label. Press firmly on the arm to cut the gel and leave in position.

Plate partitioning is recommended if only part of the plate is to be used initially or when measuring suspected high concentration samples which could (by diffusing over a wide area) result in antibody depletion occurring elsewhere on the plate. After initial use, partitioned plates should be resealed in their foil pouches and stored at 2-8°C with the gel divider(s) in place. Store partitioned plates right side up and use within four weeks.

8.3.2 Calibrator

The lyophilised calibrators should be reconstituted with the volume of the distilled water indicated on the vial labels – use the distilled water provided in the kit. Before use, all material in the bottle, including any adhering to the bung must be completely dissolved (by inversion) over a minimum period of thirty minutes. The diluted calibrators should be applied to the plates, mixing gently immediately before use. The medium and low calibrators should only be used when a calibration curve is required, as for Procedures TWO and THREE.

8.3.3 Control

The lyophilised control serum should be reconstituted with the volume of distilled water indicated on the vial label. It should be mixed gently by inversion until the contents are completely dissolved. It should then be applied to the plate(s) diluted 1/2 (i.e. 1 part in 2); for this it is recommended that 25µL of control serum is mixed with 25µL of the diluent 7% BSA.

8.3.4 Sample

Samples should be diluted 1/2 (i.e. 1 part in 2) prior to assay. To obtain adequate accuracy it is recommended that 25µL of test sample is mixed with 25µL of the diluent (7% BSA). Samples containing very high levels of C1q may require a higher dilution factor. In such cases

it is suggested that a minimum volume of 25µL of test sample is mixed with the appropriate volume of BSA. For samples having C1q concentrations below the detection limit of the plates, the following is recommended:

- i) Apply the sample undiluted
- ii) Concentrate the sample
- iii) Make a double fill of the well (see Section 8.5)

8.4 Procedures

8.4.1 Procedure ONE: RID reference table

This method does not require the construction of a calibration curve – sample concentrations corresponding to each ring diameter are read directly off the RID reference table. Rings must be allowed to develop to completion which will require a minimum diffusion time of 96 hours. The high calibrator should be run on each plate used to ensure all are performing correctly.

8.4.2 Procedure TWO: Calibration curve at completion

In this method, all three calibrators are used to produce a linear calibration curve. Rings must be allowed to develop to completion which will require a minimum diffusion time of 96 hours. To conserve wells, one calibration curve can be used for several plates of the same batch used concurrently. In such cases, the high calibrator should be run on each plate used to ensure all are performing correctly.

8.4.3 Procedure THREE: Calibration curve prior to completion

In this method, all three calibrators are used to produce a calibration curve which is non-linear, as the rings are measured before completion. The minimum recommended diffusion time is 42 hours. It is advisable that a separate calibration curve is constructed for each plate used.

8.5 Application of calibrators and samples

The calibrators, control and test samples should be mixed gently immediately before use. Fill the required number of wells with 10µL of the high calibrator using a micropipette. If Procedure TWO or THREE is being followed fill the required number of wells with the medium and low calibrators as well. The remaining wells should then be filled with 10µL of appropriately diluted test samples and controls. Plates should not be left open for long periods during calibrator test sample application, as this will cause excessive drying of the gel.

Note: For those samples suspected of containing low concentrations of C1q, a 'double fill' of the well may be made. The well is initially filled with 10µL of the sample and this is allowed to completely diffuse into the gel, which can take up to 30 minutes. The lid should be kept in place during this period. The second fill (again using 10µL) may then be made, and the plate incubated as normal. Results obtained must be corrected for the double sample volume and will be less accurate than those obtained by the normal 'single fill' procedure.

8.6 Incubation

After sample application, the lid is tightly closed and the plate stored flat at room temperature (approximately 20-24°C). It is essential that the gel is not allowed to dry out during incubation. To minimise evaporation, it is suggested that plates should either be resealed in their foil pouches or stored in a moist box (a sealed plastic box containing damp tissue paper) during incubation. The minimum incubation time for Procedure THREE is 42 hours and for complete diffusion (Procedures ONE and TWO) is 96 hours. Final ring diameters may be affected by temperature; the expected ring size for the high calibrator is 8mm (±0.3mm) when incubated at 20-24°C. Extremes of temperature should be avoided.

8.7 Quality control

The control serum should be treated exactly like a test sample, i.e. diluted 1/2. Values obtained for the control should be within ±10% of the concentration stated on the vial label.

9 RING MEASUREMENT AND RESULT PROCESSING

After the required diffusion time, ring diameters should be measured to the nearest 0.1mm, using a jewellers' eyepiece or a RID plate reader. When reading with an eyepiece, use bright side lighting and a dark background. If difficulties are experienced, view the plate macroscopically and mark the edges of the rings on the back of the plate using a needle. The distance between these marks may then be more easily measured.

Note: For Procedures ONE and TWO ring diameters must have developed to completion. If there is any doubt, rings should be remeasured after a further 24 hours to ensure there has been no increase in their diameters. The high calibrator should give a ring diameter of 8.0mm ±0.3mm at completion. If the ring diameter is outside this range, see Trouble Shooting (Section 10.3).

Procedure ONE

Read the sample concentrations directly from the RID reference table. Concentrations obtained for samples giving ring diameters greater than the high calibrator should be regarded as approximate, due to the possibility of incomplete diffusion; they may also cause local antibody depletion thereby affecting adjacent ring sizes. Such samples should preferably be diluted appropriately and retested. Samples giving ring diameters below the lower limit on the RID reference table should be retested in a more concentrated form (see Section 8.3.4). **Use of the RID reference table assumes that test samples have been applied diluted 1/2 as recommended; any change from this must be taken into account when calculating the results.**

Example:

Test sample	Dilution	Ring diameter (mm)	Table value (mg/L)	Original sample conc. (mg/L)
C1q serum A	1/2	6.4	115	115
C1q serum B	1/2	>10	>410	>410
C1q serum B (repeat)	1/4	8.0	230	460*

* Calculated as follows: Table value x Recommended Diln./Actual Diln., ie 230mg/L x (1/2)/(1/4). Note: The calibrators provided are prediluted and applied to the plate neat, not diluted 1/2. Therefore a C1q high calibrator (115mg/L) giving an 8.0mm ring is equivalent to an original sample concentration of 115 x 2 = 230mg/L, the RID reference table value.

Procedure TWO


Plot the square of the diameters of the precipitin rings formed by three calibrators versus their C1q concentrations (given on the calibrator vial label). C1q concentrations should be along the horizontal (x) axis, ring diameters squared along the vertical (y) axis. A line of best fit is drawn through the three points; the y-intercept should be in range 17-23mm². The C1q concentration is determined from the calibration curve; remember to adjust the sample concentration obtained by any dilution factor used.

Sample calculation:

C1q calibrators gave the following ring diameters on a C1q test plate at completion:

Calibrator	Conc. (mg/L)	Diameter (D) of ring (mm)	D squared (mm ²)
High	115	8.0	64.0
Medium	69	6.8	46.2
Low	11.5	4.8	23.0

A calibration curve was plotted using these results:

An unknown sample, diluted 1/2 as recommended, gave a 6.3mm diameter ring on this plate. From the above curve, this corresponds to a C1q concentration of 52mg/L. Therefore the C1q concentration in the undiluted sample = 52 x 2 = 104mg/L.

Procedure THREE


Plot the calibration curve as for procedure TWO. The graph will not be a straight line but a curve, the gradient of which decreases with increasing protein concentration. The y-intercept should be as indicated for Procedure TWO. Test sample protein concentrations are read off the calibration curve; remember to adjust the sample concentration obtained by any dilution factor used.

Sample calculation:

C1q calibrators gave the following ring diameters on a C1q plate after 42 hours:

Calibrator	Conc. (mg/L)	Diameter (D) of ring (mm)	D squared (mm ²)
High	115	7.4	54.8
Medium	69	6.8	46.2
Low	11.5	4.8	23.0

A calibration curve was plotted using these results:

An unknown sample, diluted 1/2 as recommended, gave a 6.2mm ring on this plate. From the above curve, this corresponds to a C1q concentration of 44mg/L. Therefore the C1q concentration in the undiluted sample = 44 x 2 = 88mg/L.

10 LIMITATIONS OF PROCEDURE

10.1 For Procedure ONE, results generated from ring diameters greater than the high calibrator ring diameter (i.e. 8mm) should be regarded as approximate (see Section 9). For Procedure TWO and THREE, accurate results are limited to the calibration curve between the high and low calibrator values – extrapolation beyond these points is not valid. Samples giving results outside these ranges must be diluted or concentrated as appropriate and retested (see Section 8.3.4).

10.2 FDA (USA) Information – see front page

10.3 TROUBLE SHOOTING

Problem	Possible causes(s)	Suggested action(s)
A. No ring for:		
1. Calibrator(s)	Calibrator omitted.	Repeat assay.
2. Test sample	i) Sample omitted.	Repeat assay.
	ii) Concentration too high/low.	Dilute/concentrate and reassay.
3. Calibrator(s) and test samples	Plate deterioration	a) Storage damage. Repeat assay using new plate. b) Product expired. Repeat assay using new plate/kit.

Problem	Possible causes(s)	Suggested action(s)
B. Oversize rings for:		
1. High calibrator (more than 8.3mm)	i) Inaccurate ring measurement.	Remeasure using eyepiece or a RID plate reader.
	ii) Incorrect volume applied.	Check 10µL volume applied.
	iii) Inaccurate volume applied.	a) Micropipette malfunction – check operation and repeat assay. b) Poor technique – repeat assay.
	iv) Inaccurate calibrator reconstitution.	a) Pipette malfunction – check operation and calibration, then repeat assay using new calibrator. b) Poor technique – repeat assay using new calibrator.
	v) Partial evaporation of reconstituted calibrator on storage.	Repeat assay using new calibrator/kit.
	vi) Plate deterioration.	a) Storage damage. Repeat assay using new plate. b) Product expired. Repeat assay using new kit.
	vii) Local antibody depletion due to adjacent high concentration test samples.	Dilute the sample(s) responsible and repeat assay using new plate.
	viii) Incubation temperature too high (see Section 8.6).	Repeat assay, incubating at 20-24°C.
2. Test samples (above acceptable range - see section 10.1).	i) Concentration too high.	Dilute and reassay.
	ii) Incorrect volumes applied.	Check 10µL volume applied.
C. Undersized rings for:		
1. High calibrator (less than 7.7mm)	i) Inaccurate ring measurement.	As for B1 above
	ii) Incorrect volume applied.	
	iii) Inaccurate volume applied.	
	iv) Inaccurate calibrator reconstitution	
	v) Calibrator deterioration.	a) Storage damage. Repeat assay using new calibrator. b) Product expired. Repeat assay using new kit.
	vi) Incubation temperature too low (see Section 8.6).	Repeat assay, incubating at 20-24°C.
2. Test samples (below acceptable range – see Section 10.1).	i) Concentration too low.	See section 8.3.4 and repeat assay.
	ii) Incorrect volume applied.	Check 10µL volume applied.
D. Double/multiple rings		
	i) Non-specific precipitation close to well (due to PEG in gel).	Read outer ring.
	ii) Poor sample application.	Repeat assay.
	iii) Calibrator deterioration.	a) Storage damage. Repeat assay using new calibrator. b) Product expired. Repeat assay using new kit.
	iv) Sample deterioration.	Reassay using fresh sample.
E. Non-circular rings		
	i) Poor sample application.	Repeat assay.
	ii) Gel dried out before use.	a) Storage damage. Repeat assay using new plate. b) Product expired. Repeat assay using new plate/kit.
	iii) Gel dried out during sample application or incubation.	Repeat assay minimising the time the plate is left open. Incubate with lid on tight in a moist box or sealed foil pouch.
	iv) Local antibody depletion (due to high concentration samples on the plate).	Dilute samples and repeat assay.
F. Cloudy gel		
	i) Plate has been frozen.	Repeat assay using new plates. Review storage.
	ii) Gel dried out before use.	As for E(ii) above.
	iii) Gel dried out during sample application or incubation.	As for E(iii) above.
G. Weak, pitted gel		
	Plate has been frozen.	Repeat using new plate. Review storage.
H Poor calibration curve		
1. Curve non-linear (Procedure TWO)	i) Incomplete diffusion.	Incubate for further 24 hours and remeasure the rings.
	ii) Calibrator rings under/oversize.	As for B1 or C1 above. (Similar explanations apply to the medium and low calibrators).
	iii) Calibration curve constructed incorrectly.	Check calibration curve construction.
2. y-intercept out-of-range (Section 9)	i) Calibrator rings under/oversize.	As for B1 or C1 above. (Similar explanations apply to the medium and low calibrators).
	ii) Calibration curve constructed incorrectly.	Check calibration curve construction.

10.4 Diagnosis cannot be made and treatment must not be initiated on the basis of C1q measurements alone. Clinical history and other laboratory findings must be taken into account.

10.5 If an unexpected result is obtained, the assay should be repeated, preferably with a fresh sample.

If a problem cannot be resolved, please refer to supplier.

11 EXPECTED VALUES

The following results were obtained using this kit using individual blood donors.

	No of samples	Mean (mg/L)	Median (mg/L)	Standard deviation	95 percentile Range
Normal males	62	163	155	35	118-238
Normal females	60	158	151	33	118-244
Active SLE	14	117	122	52	33-209
Inactive SLE	17	128	126	28	93-183

The data provided above has been obtained from limited numbers of British blood donors and clinical patients and is intended for guidance purposes only. It is strongly recommended that each user should generate his/her own C1q concentration ranges for appropriate clinical conditions.

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

The precision (repeatability) of this kit is expressed as the mean and the percentage coefficient of variation (CV), which had been determined using human serum preparation containing high, medium and low concentrations of C1q. All analyses were performed in our laboratory. Each value was calculated from 10 measurements (duplicate determinations on five separate plates from a typical batch) unless otherwise stated. For Procedures ONE and TWO, rings were measured after 96 hours. For Procedure THREE, rings were read after 42 hours.

Sample pool C1q	Procedure 1		Procedure 2		Procedure 3	
	Mean conc. mg/L	CV	Mean conc. mg/L	CV	Mean conc. mg/L	CV
High	191.6	1.76%	189.0	2.23%	178.2	7.62%
Medium	125.7	4.77%	119.2	4.63%	118.0	7.85%
Low	44.96	6.90%	39.59	7.82%	42.55	8.55%

12.2 Within-plate and inter-batch variation:

The within-plate variation is expressed as the mean ± standard deviation of determinations of CV made using 3 plates from separate batches. Six measurements were made per plate, using a human serum pool as the sample.

The inter-batch variation is expressed as the CV of mean diameter values obtained from recent batches of plates. The mean diameter for each batch was calculated using the ring diameter at completion obtained using a serum pool as the sample, applied to two plates from each batch (six ring measurements per plate).

C1q	Within-plate variation	Interbatch variation
	Mean CV% ± SD	CV (%)
	0.79 ± 0.16 (N=3)	0.23 (N=3)

15 RID REFERENCE TABLE

RID reference table for human C1q Concentrations in mg/L

Diameter of ring (mm)	Conc.
4.5	11.0
4.6	15.5
4.7	20.2
4.8	25.0
4.9	29.8
5.0	34.8
5.1	39.8
5.2	45.0
5.3	50.2
5.4	55.6
5.5	61.1
5.6	66.6
5.7	72.3
5.8	78.0
5.9	83.9
6.0	89.8
6.1	95.9
6.2	102
6.3	108
6.4	115
6.5	121
6.6	128
6.7	134
6.8	141
6.9	148
7.0	155
7.1	162
7.2	169
7.3	176
7.4	184
7.5	191
7.6	199
7.7	206
7.8	214
7.9	222
8.0	230
8.1	238
8.2	246
8.3	254
8.4	263
8.5	271
8.6	280
8.7	289
8.8	297
8.9	306
9.0	315
9.1	324
9.2	333
9.3	343
9.4	352
9.5	361
9.6	371
9.7	381
9.8	390
9.9	400
10.0	410

Note: The above values assume that test samples are applied diluted 1/2 in 10µL volumes. The high calibrator should give a ring diameter of 8.0 ± 0.3mm at completion when incubated at 20-24°C.

13 BIBLIOGRAPHY

13.1 Cooper, NR, (1985). The classical complement pathway: Activation and regulation of the first complement component. *Adv. Immunol.* **37**, 151-216

13.2 Aritaud, GJ *et al* (1987). A functional model of the human C1 complex. *Immunol. Today*, **8**, 106-111

13.3 Ross, SC & Densen, P (1984). Complement deficiency states and infection: Epidemiology, pathogenesis and consequences of neisserial and other infections in an immune deficiency. *Medicine*, **63**, 243-273

13.4 Schifferli, J A *et al* (1986). The role of complement and its receptor in the elimination of immune complexes. *New Eng. J.* **315**, 488-495

13.5 Fahey, JL & McKelvey, EM (1965). Quantitative determination of serum immunoglobulins in antibody-agar plates. *J. Immunol.*, **94**, 84-90.

13.6 Mancini, G, Vaerman, J P *et al* (1964). Protides of the biological fluids (XI Colloquium). Peters H. (ed), Amsterdam, Elsevier Publishing Co., p 370.

13.7 Mancini, G, Carbonara, A O *et al* (1965). Immunochemical quantification of antigens by single radial immunodiffusion. *Immunochem.* **2**, 235-254.

14 SUMMARY OF PROCEDURE

- 14.1 Select Procedure ONE, TWO or THREE. Procedure THREE must be used if results are required quickly.
- 14.2 Reconstitute calibrator(s) and control with the distilled water provided.
- 14.3 Prepare 1/2 sample and control dilutions with diluent (7% BSA) provided.
- 14.4 Allow condensation to evaporate from RID plate(s).
- 14.5 Apply calibrator(s), control and samples to RID plate(s) in 10µL volumes.
- 14.6 Replace lid and incubate at room temperature (approximately 20-24°C) for fixed time period (minimum 42 hours) (Procedure THREE) or until rings are complete (minimum 96 hours) (Procedure ONE and TWO).
- 14.7 Measure the ring diameters.
- 14.8 Read results off RID Reference Table (Procedure ONE) or plot calibration curve and read off results (Procedures TWO and THREE).

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13 BIBLIOGRAPHIE

- Cooper, NR, (1985). The classical complement pathway: Activation and regulation of the first complement component. Adv. Immunol. **37**, 151-216
- Arlaud, GJ *et al* (1987). A functional model of the human C1 complex. Immunol. Today, **8**, 106-111
- Ross, SC & Densen, P (1984). Complement deficiency states and infection: Epidemiology, pathogenesis and consequences of neisserial and other infections in an immune deficiency. Medicine, **63**, 243-273
- Schifferli, J A *et al* (1986). The role of complement and its receptor in the elimination of immune complexes. New Eng. J. **315**, 488-495
- Fahey, JL & McKelvey, EM (1965). Quantitative determination of serum immunoglobulins in antibody-agar plates. J. Immunol., **94**, 84-90.
- Mancini, G, Vaerman, J P *et al* (1964). Protides of the biological fluids (XI Colloquium). Peters H. (ed), Amsterdam, Elsevier Publishing Co., p 370

- 13.7 Mancini, G, Carbonara, A O *et al* (1965). Immunochemical quantification of antigens by single radial immunodiffusion. Immunochem. **2**, 235-254

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13 BIBLIOGRAPHIE

- Cooper, NR, (1985). The classical complement pathway: Activation and regulation of the first complement component. *Adv. Immunol.* **37**, 151-216
- Arlaud, GJ *et al* (1987). A functional model of the human C1 complex. *Immunol. Today*, **8**, 106-111
- Ross, SC & Densen, P (1984). Complement deficiency states and infection: Epidemiology, pathogenesis and consequences of neisserial and other infections in an immune deficiency. *Medicine*, **63**, 243-273
- Schifferli, J A *et al*. (1986). The role of complement and its receptor in the elimination of immune complexes. *New Eng. J.* **315**, 488-495
- Fahey, JL & McKelvey, EM (1965). Quantitative determination of serum immunoglobulins in antibody-agar plates. *J. Immunol.*, **94**, 84-90.
- Mancini, G, Vaerman, J P *et al*. (1964). Protides of the biological fluids (XI Colloquium). Peters H. (ed), Amsterdam, Elsevier Publishing Co., p 370
- Mancini, G, Carbonara, A O *et al* (1965). Immunochemical quantification of antigens by single radial immunodiffusion. *Immunochem.* **2**, 235-254

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13 BIBLIOGRAFIA

- Cooper. NR. (1985). The classical complement pathway: Activation and regulation of the first complement component. Adv. Immunol. **37**, 151-216
- Arlaud, GJ *et al* (1987). A functional model of the human C1 complex. Immunol. Today, **8**, 106-111
- Ross, SC & Densen, P (1984). Complement deficiency states and infection: Epidemiology, pathogenesis and consequences of neisserial and other infections in an immune deficiency. Medicine, **63**, 243-273
- Schifferli, J A *et al*. (1986). The role of complement and its receptor in the elimination of immune complexes. New Eng. J. **315**, 488-495
- Fahey, JL & McKelvey, EM (1965). Quantitative determination of serum immunoglobulins in antibody-agar plates. J. Immunol., **94**, 84-90.
- Mancini, G, Vaerman, J P *et al*. (1964). Protides of the biological fluids (XI Colloquium). Peters H. (ed), Amsterdam, Elsevier Publishing Co., p 370.
- Mancini, G, Carbonara, A O *et al* (1965). Immunochemical quantification of antigens by single radial immunodiffusion. Immunochem. **2**, 235-254.

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